

Instructions to setup Iriun Webcam (Windows)

1. Download Iriun Webcam app on your android phone from playstore.
2. Download the iriun webcam app client for windows from the given source in the installation page.
3. Go to your smartphone settings and enable USB debugging follow the below link for details:
 - a. <https://developer.android.com/studio/debug/dev-options>
4. Connect your smartphone to your computer via USB.
5. Open the Iriun webcam app on your smartphone and the Iriun webcam client on your computer.
6. Allow the client to display the camera feed from your smartphone.
7. You should now see the video feed from your Android smartphone on your PC.
8. Adjust the positioning of your camera as instructed in the setting up tutorial
9. **Keep this video feed tab minimized and running all the time while you perform exercises using the Mhandrehab app.**
10. In case the feed is not properly setup you will see an animated cartoon on a black screen notifying you that video feed is unavailable.
11. Please reach out to program supervisors in case any issues.